

SHADE COMBINATION CHART

		A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4	A0	B0	
FOUNDATION OPAQUE																				
	MARGIN OPAQUE																			
OPAGUE	OPAQUE ORIGINAL SHADE																			
	OPAQUE CLASSICAL SHADE	OA1-C	OA2-C	OA3-C	OA3.5-C	OA4-C														
	OPAQUE DENTIN, INTENSIVE*1 ORIGINAL SHADE	ODA2	ODA3	ODA3.5	ODA4	ODI3								ODI5						
	OPAQUE DENTIN, INTENSIVE*1 CLASSICAL SHADE	ODA2-C	ODA3-C	ODA3.5-C	ODA4-C															
DENTIN	OPAQUE DENTIN ORIGINAL SHADE	ODA1	ODA2	ODA3	ODA3.5	ODA4														
	OPAQUE DENTIN CLASSICAL SHADE	ODA1-C	ODA2-C	ODA3-C	ODA3.5-C	ODA4-C														
	DENTIN ORIGINAL SHADE	DA1	DA2	DA3	DA3.5	DA4														
	DENTIN CLASSICAL SHADE	DA1-C	DA2-C	DA3-C	DA3.5-C	DA4-C														
ENAMEL	ENAMEL																			
	ENAMEL INTENSIVE																			
TRANSLUCENT	CERVICAL TRANSLUCENT	*2				*2			*2											

*1: OPAQUE DENTIN and OPAQUE DENTIN INTENSIVE are used as cervical colors.
 *2: The use of ENAMEL INTENSIVE (E11) will give almost the same effects as CERVICAL TRANSLUCENT.