

GRADA CHARACTERIZATION

SHOULDER DENTIN	SD2	SD3	SD4	SD5	SD7	SD8										
	OD11	OD12	OD13	OD14	OD15	OD16										
	ENAMEL INTENSIVE	E10	E11	E13	E15											
	HALO ENAMEL	HE1														
	PEARL ENAMEL	PE1	PE3													
		T0	T1	T2	T4	T5										
	CERVICAL TRANSLUCENT	CT2	CT3	CT4												
		MAMELON STAIN	MS2	MS3	MS5											
INTENSIVE COLOR		IC0	IC1	IC2	IC3	IC4	IC5	IC6	IC7	IC8	IC9	IC10	IC11	IC12	IC13	IC14
		FLOWABLE COMPOSITE	F-WT	F-WE												

IRRADIATION TIME CHART

	LABOLIGHT LV-III. II		STEPLIGHT SL-I
	Pre-cure	Final Cure	
FOUNDATION OPAQUE			
OPAQUE			
MARGIN OPAQUE	1 min.	—	—
GUM OPAQUE			
GUM OPAQUE MODIFIER			
DENTIN			
OPAQUE DENTIN			
OPAQUE DENTIN INTENSIVE			
SHOULDER DENTIN			
ENAMEL			
ENAMEL INTENSIVE			
HALO ENAMEL	30 sec.	3 min.	10 sec.*
PEARL ENAMEL			
TRANSLUCENT			
CERVICAL TRANSLUCENT			
GUM			
GUM MODIFIER			
GUM TRANSLUCENT			
INTENSIVE COLOR			
MAMELON STAIN	30 sec.	—	10 sec.*
FLOWABLE COMPOSITE			

*For one surface of a single crown